

	TECHNISCHES DATENBLATT	<b>Code:</b> TDS_095
	KÄLTESPRAY PROCLINIC	<b>Datum:</b> 26. Juni 2019
	PROCLINIC	<b>Version:</b> 01

<b>Produkt</b>	<b>KÄLTESPRAY PROCLINIC</b>	
<b>Artikel-Nr.</b>	01924 01923	
<b>Hersteller</b>	PRISMAN GmbH. Otto-Hahn-Ring 6-18 DE-64653 Lorsch Tel: +49 6251 866980-0	
<b>Klassifizierung</b>	Medizinprodukt Klasse I, gemäß der Richtlinie 93/42/EWG über Medizinprodukte.	
<b>Produktbeschreibung</b>	In zwei angenehmen Duftrichtungen erhältlich, Orange und Minze. Mit Hilfe des speziellen Applikators lässt es sich präzise applizieren. Bei jeder Applikation wird damit auch Spray eingespart.  Abkühlung auf bis zu -40 °C.	
<b>Technische Daten</b>	01924: Minzeduft  01923: Orangenduft	
<b>Packungen</b>	Dose mit 200 ml	
<b>Haltbarkeit</b>	3 Jahre	
<b>Bestimmungsgemäßer Gebrauch</b>	Kältespray zur Vitalitätsprüfung der Pulpa und Abkühlung von Pellets, Glasanmischplatten, etc.	
<b>Sicherheitshinweise</b>	Von Wärme, heißen Oberflächen, Funken, Flammen und anderen Zündquellen fernhalten. In der Umgebung nicht rauchen. Auch nach Gebrauch nicht durchbohren oder verbrennen. Nicht in offene Flammen oder andere Zündquellen sprühen. Vor Sonnenbestrahlung und Temperaturen über 50 ° C schützen.	
<b>Transport und Lagerung</b>	Kühl und trocken aufbewahren. Behälter dicht verschlossen lassen. Vor Wärme und direkter Sonnenbestrahlung schützen. Senkrecht stehend lagern.	

	TECHNICAL DATA SHEET	CODE: TDS_095
	REFRIGERATING SPRAY PROCLINIC	DATE: 26/JUN/19
	PROCLINIC	VERSION: 01

<b>PRODUCT</b>	REFRIGERATING SPRAY PROCLINIC	
<b>REFERENCE</b>	01924 01923	
<b>COMPANY</b>	PRISMAN GmbH. Otto-Hahn-Ring 6-18 64653 Lorsch Tel: +49 6251 866980-0	
<b>CLASSIFICATION</b>	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
<b>PRODUCT DESCRIPTION</b>	Available in two pleasant fragrances such as orange and mint, it is applied with precision thanks to its special applicator that also saves spray in each application.  Cool up to - 40 ° C.	
<b>SPECIFICATIONS</b>	01924: Mint Fragrance  01923: Fragrance Orange	
<b>PACKAGING</b>	Bottle of 200 mL.	
<b>SHELF LIFE</b>	3 years	
<b>INTENDED OF USE</b>	Coolant spray for pulp vitality test and cooling of pellets, glass tiles, etc.	
<b>SECURITY</b>	Keep away from heat, hot surfaces, sparks, flames and other sources of ignition. Do not smoke nearby. Do not pierce or burn, even after use. Do not spray on an open flame or other source of ignition. Protect from sunlight. Do not expose to temperatures above 50 ° C.	
<b>TRANSPORT AND STORAGE CONDITIONS</b>	Store in a cool place. Keep container tightly closed. Store in a cool, dry place in tightly closed containers. Protect from heat and direct sunlight. Store in vertical position.	